

Neoss Clinical Organisers


Simplicity and ease of use are two of the most important parts of using a dental implant system. It may be said that the instrument tray is the key to this being used for organisation, storage, sterilisation and information.

The Neoss Clinical Organisers are designed as three interlocking parts ('a jigsaw') for surgery, instruments and layout. Made of highly durable silicone it is easily cleaned and sterilised.

The surgical part of the organiser (one for Neoss ProActive® implants and one for Neoss Tapered implants) offers clear markings for drill selection and depth on one side and offers storage for instruments and drills during sterilisation on the other. The mid section may be used in combination with the other parts or alone for prosthodontics.

The layout section on the left provides wells for implant storage, cover and abutment screws on one side and prosthetic components crowns and bridges on the other.

The Neoss Clinical Organisers – Simplicity Reinvented



Neoss Clinical Organiser features:

- A unique three part storage and instrument tray solution which can be used in combination or individually.
- Extremely competitive costing enables use of multiple organisers for serial cases.
- Clear identification of drill lengths and selection, with horizontal instrument storage for sterilisation.
- Instrument section usable alone or in combination.
- Storage section for easy implant selection and convenient arrangement for prosthetic component storage.
- Unrivalled ease of use and cleaning – 'The Nurses Friend'.


Assortment

Item no.	Description
51150	Neoss Clinical Organiser
51151	Neoss Clinical Organiser, Tapered

Note: Neoss system tray, 51131, includes customer choice of 51150 or 51151

Overview

Surgical setup



Prosthetic setup



Storage and sterilisation

